



**DONGSHENG BIOTECH**

## M-MLV Reverse Transcriptase

**#R1042**

10,000U

|                       |                       |             |
|-----------------------|-----------------------|-------------|
| <b>Concentration:</b> | 200U/ $\mu$ l         |             |
| <b>Contents:</b>      | M-MLV                 | 50 $\mu$ l  |
|                       | 5xfirst-strand buffer | 200 $\mu$ l |

**Store at -20°C**

For research use only

In total 2 vials.

### ***Description***

Moloney Murine Leukemia Virus Reverse Transcriptase (M-MLV RT) uses single-stranded RNA or DNA in the presence of a primer to synthesize a complementary DNA strand. This enzyme is isolated from *E. coli* expressing a portion of the pol gene of M-MLV on a plasmid. The enzyme is used to synthesize first-strand cDNA up to 5 kb.

### ***Unit Definition***

One unit incorporates 1 nmole of dTTP into acid-precipitable material in 10 min at 37 ° C using poly(A) • oligo(dT) 25 as template-primer.

### ***Source***

Purified from an *E. coli* strain expressing a recombinant clone.

### ***Storage Buffer***

20mM Tris-HCl (pH 7.5), 200mM NaCl, 0.1mM EDTA, 1mM DTT, 0.01% NP-40, 50% glycerol

### **PRODUCT USE LIMITATION.**

This product is developed, designed and sold exclusively *for research purposes and in vitro use only*. The product was not tested for use in diagnostics or for drug development, nor is it suitable for administration to humans or animals.

## ***5xfirst-Strand Buffer***

250mM Tris-HCl (pH 8.3 at 25°C), 375mM KCl, 15mM MgCl<sub>2</sub>  
50mM DTT

## ***Applications***

Generation of first strand cDNA for use in:

- PCR, see Protocol for First-strand cDNA Synthesis;
- Second strand cDNA synthesis.
- DNA labeling.
- Real-time PCR;
- Analysis of RNA by primer extension.

## ***Protocol***

### **I. First-Strand cDNA Synthesis Using M-MLV RT**

A 20- $\mu$ l reaction volume can be used for 1ng–5 $\mu$ g of total RNA or 1–500ng of mRNA.

**1. Add the following reagents into a sterile, nuclease-free tube on ice in the indicated order:**

|                    |  |                           |
|--------------------|--|---------------------------|
| Template RNA       | poly(A) mRNA<br>or specific RNA  | 1 to 500ng<br>1-5 $\mu$ g |
| Prime              | oligo (dT) <sub>15</sub> primer(50 $\mu$ M)<br>or Random hexamer<br>primer(50 $\mu$ M) | 1 $\mu$ l<br>1 $\mu$ l    |
| DEPC-treated water |  | to 13.4 $\mu$ l           |
| Total volume       |  | 13.4 $\mu$ l              |

**2. Mix gently, centrifuge briefly and incubate at 70°C for 5 min. Chill on ice, spin down and place the vial back on ice.**

**3. Prepare the following cDNA Synthesis Mix, add the following components in the indicated order:**

|                        |             |
|------------------------|-------------|
| 5x first-strand buffer | 4 $\mu$ l   |
| dNTPs(10 mM each)      | 1 $\mu$ l   |
| RNasin (40U/ $\mu$ l)  | 0.6 $\mu$ l |
| M-MLV                  | 1 $\mu$ l   |

**4. Mix gently and centrifuge**

**5. For oligo(dT)<sub>15</sub>, incubate for 60 min at 42°C. For random hexamer primed synthesis, incubate for 60 min at 37°C .**

**6. Terminate the reaction by heating at 70°C for 5 min.**

The reverse transcription reaction product can be used immediately in second strand cDNA synthesis reactions or stored at -20 ° C for less than a week. For longer storage, -70 ° C is recommended.

## **II. PCR Reaction**

The product of the first strand cDNA synthesis can be used directly in PCR or qPCR. The volume of first strand cDNA synthesis reaction mixture should not comprise more than 1/10 of the total PCR reaction volume. Normally, 2  $\mu$ l of the first strand cDNA synthesis reaction mixture is used as template for subsequent PCR in 50  $\mu$ l total volume.

## **Quality Control**

This product has passed the following quality control assays: SDS polyacrylamide gel analysis for purity; functional absence of endodeoxyribonuclease, 3' and 5' exodeoxyribonuclease, and ribonuclease activities, yield and length of cDNA product.